

Link do produktu: <https://higienadent.com/alginat-do-wyciskow-ortodontycznych-lascod-alginor-450g-p-1911.html>



Alginate for dental impressions Lascod Alginor / 450g

Cena	17,00 zł
Dostępność	Niedostępny
Numer katalogowy	KAT00212
Producent	Lascod

Opis produktu

Lascod Alginor alginate for dental impressions with chromatic phase indicator : violet -> pink -> green. Its high elasticity and resistance to tearing allow for perfect execution of impressions in the presence of visible undercuts and already installed brackets.

Features of the product:

- High elasticity
- Alginate class A type 1
- Easy to mix
- Thixotropic
- Very fast setting time in the oral cavity 30"
- High accuracy: 20 µm
- Additional dimensional stability (above 168 hours)
- Compatible with alginate disinfectant for impressions
- Minty aroma

Properties: chromatic alginate for dental impressions with high elasticity and resistance to tearing for accurate impressions with deep undercuts and brackets. It can be used for joining impressions for:

- master models
- research models
- standard dental models

Packaging: 450g